

Media Contact: Tom Greenwell **Tel:** +44 (0)870 991 4044 **Email:** tom.greenwell@analytik.co.uk

Analytik launch new FieldSpec HandHeld 2™ portable spectroradiometer

Cambridge, UK – Analytik are delighted to announce the official UK and Ireland launch of the new FieldSpec HandHeld 2™ portable spectroradiometer from ASD, for remote sensing research, oceanography and environmental science.

The HandHeld 2™ is the latest in compact remote sensing technology from ASD using cutting-edge, near-infrared (NIR) spectroscopy to instantly and non-destructively make material measurements. Self-contained data collection with significant spectrum storage capacity and an on-board colour display facilitate more efficient field work by eliminating the need to be continuously connected to a computer.

Other enhancements include ergonomically designed controls, laser targeting, standard AA battery power and a rugged enclosure to ensure durability and easy transport.

“Building on the company’s excellence in portable NIR, the HandHeld 2 brings fast, laboratory-grade spectroscopy to field research,” said ASD CEO David Rzasa. “Well-thought out enhancements, like laser targeting and onboard data collection, make this instrument a top choice for field researchers and scientists across the world.”

The HandHeld 2™ provides the performance required for a broad range of applications with a wavelength range of 325-1075 nm, accuracy of +/- 1.0 nm, 3 nm resolution, standard 25° field-of-view, proprietary DriftLock™ compensation for dark current and fast scan times. The unit weighs less than three pounds, measures 215x140x90 mm, and uses four rechargeable or standard AA batteries. The HandHeld 2™ comes standard with a Pelican case that meets FAA requirements for carry-on baggage.



Attachment: The new FieldSpec HandHeld 2 field spectroradiometer

To find out more about the FieldSpec HandHeld 2™, visit <http://www.asdi.com/products/handheld-2-portable-spectroradiometer>

About Analytik:

Analytik are specialised consultants and suppliers of unique, high quality analytical instrumentation throughout UK and Ireland. Innovative solutions from leading global instrumentation manufacturers, coupled with unparalleled customer support have enabled Analytik to forge a reputation for excellence since forming in 2003.

###

For further information contact:

analytikLtd

Barn B, 2 Cygnus Business Park, Middle Watch, Swavesey, Cambridge, CB24 4AA

T: +44 (0)870 991 4044 **F:** +44 (0)870 135 2488 **E:** info@analytik.co.uk **www.analytik.co.uk**